

Data Sheet

Model: TMH15(Z3)/50/450

1500 °C Maximum Temperature,
Three Zone Horizontal Tube Furnace



INTRODUCTION

The TMH15(Z3) furnace is a bench mounted tube furnace ideal for most general laboratory thermal processing applications.

The TMH15(Z3)/50 Three Zone Tube furnace has an internal tube diameter of 50 mm. TMH15(Z3)/50 range offers a choice of two heated zone lengths of 450mm & 610mm.

TMH15(Z3)/50/450 comes with a heated zone length of 450 mm.



TMH15(Z3)/50/450

SPECIFICATIONS

Maximum Temperature: 1500 °C

Maximum Continuous Temperature: 1500 °C

- | The TMH(Z3) furnaces are designed to give a longer uniform centre zone temperature than that of the single zone tube furnace models.
- | All TMH(Z3) models are controlled by retransmission of set point from the centre zone controller to the end zone controllers.
- | This system provides a longer uniform zone temperature than that achieved by the use of single zone furnace of the same length.
- | Protective outer mesh covers for improved operator safety.
- | Top grade silicon carbide rods in all 3 heated zones.
- | Energy efficient, high quality, low thermal mass insulation.
- | A work tube is not supplied as an integral part of the furnace and therefore needs to be ordered with the furnace as it is an essential accessory.
- | High end Microprocessor PID controller.
- | Temperature Sensor: "R" type Thermocouple.
- | Horizontal models are supplied as standard with controls in the base.

External Dimensions (mm): H x W x D **Net Wt.**
(Nominal) TMH15(Z3)/50/450 715 x 1000 x 460 55 kg

- | Supply: 400V– 3 Phase
- | Energy rating: 4.5 kW (Average)
22.0 kW (Peak)

OPTIONAL FEATURES

- | Work tubes of various materials, lengths and diameters for use in the furnace.
- | The work tubes are available for containment of atmosphere or protection against process contaminants.
- | A variety of triple flange gas tight end seals for work tubes to allow processing under Vacuum or gas atmospheres.
- | Over Temperature protection controller.
- | Multi segment, multi program storage controllers.

Vertical versions are available at extra cost, which will have separate temperature console

Elite Thermal Systems Limited

Elite Court, 6 Stuart Road, Market Harborough, Leicestershire LE16 9PQ, UK
Tel: +44 (0)1858 469834 | E-mail: contact@elitefurnaces.com | Website: www.elitefurnaces.com