

Data Sheet

Model: TMH17(Z3)/50/450

1700 °C Maximum Temperature,
Three Zone Horizontal Tube Furnace



TMH17(Z3)/50/450

INTRODUCTION

The TMH17(Z3) furnace is a floor mounted tube furnace ideal for most general laboratory thermal processing applications.

The TMH17(Z3)/50 Three Zone Tube furnace has an internal tube diameter of 50 mm. TMH17(Z3)/50 range offers a choice of two heated zone lengths of 450mm & 610mm.

TMH17(Z3)/50/450 comes with a heated zone length of 450 mm.

SPECIFICATIONS

Maximum Temperature: 1700 °C

Maximum Continuous Temperature: 1700 °C

- | The TMH(Z3) furnaces are designed to give a longer uniform centre zone temperature than that of the single zone tube furnace models of the same length.
- | All TMH(Z3) models are controlled by retransmission of set point from the centre zone controller to the end zone controllers.
- | Molybdenum Disilicide single loop elements mounted along the front and rear walls of the furnace chamber.
- | Energy efficient, high quality, low thermal mass insulation.
- | A work tube is not supplied as an integral part of the furnace and therefore needs to be ordered with the furnace as it is an essential accessory.
- | Temperature Sensor: 'B' type thermocouple.
- | High end Microprocessor PID controller.
- | Over temperature protection is included in the standard specification.
- | Horizontal models are supplied as standard with controls in the base.

External Dimensions (mm): H x W x D **Net Wt.**
(Nominal) TMH17(Z3)/50/450 1550 x 835 x 650 280 kg

| Supply: 230V– Single Phase

| Energy rating: 9.0 kW

PHASE CONTROL

Phase angle fired thyristor units operating in conjunction with low voltage secondary isolating transformers providing the correct parameters for the heating elements.

OPTIONAL FEATURES

- | Work tubes of various materials, lengths and diameters for use in the furnace.
- | The work tubes are available for containment of atmosphere or protection against process contaminants.
- | A variety of triple flange gas tight end seals for work tubes to allow processing under gas atmospheres.
- | Multi segment, multi program storage controller.

Vertical versions are available at extra cost, which will have separate temperature console

Elite Thermal Systems Limited

Elite Court, 6 Stuart Road, Market Harborough, Leicestershire LE16 9PQ, UK

Tel: +44 (0)1858 469834 | E-mail: contact@elitefurnaces.com | Website: www.elitefurnaces.com