

Data Sheet

Model: TMV15(Z3)/50/610

1500 °C Maximum Temperature,
Three Zone Vertical Tube Furnace



INTRODUCTION

The TMV15 furnace is a vertical-mounted tube furnace with a separate control console, ideal for most general laboratory thermal processing applications.

The TMV15(Z3)/50 Three Zone Tube furnace has an internal tube diameter of 50 mm. TMV15(Z3)/50/610 comes with a heated zone length of 610 mm.

SPECIFICATIONS

Maximum Temperature: 1500 °C

Maximum Continuous Temperature: 1500 °C

- | The TMV(Z3) furnaces are designed to give a longer uniform centre zone temperature than that of the single zone tube furnace models.
- | All TMV models are controlled by retransmission of set point from the centre zone controller to the end zone controllers.
- | This system provides a longer uniform zone temperature than that achieved by the use of single zone furnace of the same length.
- | Protective outer mesh covers for improved operator safety.
- | Top grade silicon carbide rods in all 3 heated zones.
- | Energy efficient, high quality, low thermal mass insulation.
- | Temperature Sensor: Thermocouple, type "R".
- | A work tube is not supplied as an integral part of the furnace and therefore needs to be ordered with the furnace as it is an essential accessory.
- | High end Microprocessor PID controller.
- | Horizontal models are supplied as standard with controls in the base.



TMV15(Z3)/50/610

| **External Dimensions (mm):** H x W x D **Net Wt.**
(Nominal) TMV15(Z3)/50/610

| Supply: 400V– 3 Phase

| Energy rating: 5.0 kW (Average)
23.0 kW (Peak)

OPTIONAL FEATURES

- | Work tubes of various materials, lengths and diameters for use in the furnace.
- | The work tubes are available for containment of atmosphere or protection against process contaminants.
- | A variety of end seals for work tubes to allow processing under Vacuum or gas atmospheres.
- | Over Temperature protection controller.
- | Multi segment, multi program storage controllers.

Elite Thermal Systems Limited

Elite Court, 6 Stuart Road, Market Harborough, Leicestershire LE16 9PQ, UK
Tel: +44 (0)1858 469834 | E-mail: contact@elitefurnaces.com | Website: www.elitefurnaces.com