

Data Sheet

Model: TSH18/75/300

1800 °C Maximum Temperature,
Single Zone Horizontal Tube Furnace



INTRODUCTION

The TSH18 furnace is a floor mounted tube furnace ideal for most general laboratory thermal processing applications.

The TSH18/75 Single Zone Tube furnace can accept a work tube with internal tube diameter of 75 mm. TSH18/75 range offers a choice of three heated zone lengths of 300mm, 450mm and 600mm.

TSH18/75/300 comes with a heated zone length of 300 mm.

SPECIFICATIONS

Maximum Temperature: 1800 °C

Maximum Continuous Temperature: 1800 °C

- | These models are heated on both sides of the chamber by molybdenum Tungsten Disilicide elements.
- | This furnace design requires the use of a separate work tube of a grade suitable for the maximum temperature rating of 1800 °C.
- | A work tube is not supplied as an integral part of the furnace and therefore needs to be ordered with the furnace as it is an essential accessory.
- | Temperature Sensor: 'B' type Thermocouple.
- | Energy efficient, high quality, low thermal mass insulation.
- | High end Microprocessor PID controller.
- | Over temperature protection controller.
- | Horizontal models are supplied as standard with controls in the base.



TSH18/75/300

External Dimensions (mm): H x W x D **Net Wt.**
(Nominal) TSH18/75/300 1500 x 675 x 650 220 kg

| Supply: 230V– 1 Phase – 50/60Hz

| Energy rating: 6.0 kW

PHASE CONTROL

Phase angle fired thyristor units operating in conjunction with low voltage secondary isolating transformers provide the correct parameters for the heating elements.

OPTIONAL FEATURES

- | Work tubes of various materials, lengths and diameters for use in the furnace.
- | The work tubes are available for containment of atmosphere or protection against process contaminants.
- | A variety of triple flange gas tight end seals for work tubes to allow processing under gas atmospheres.
- | Multi segment, multi program storage controllers.
- | Thermal Radiation Plugs.

Elite Thermal Systems Limited

Elite Court, 6 Stuart Road, Market Harborough, Leicestershire LE16 9PQ, UK
Tel: +44 (0)1858 469834 | E-mail: contact@elitefurnaces.com | Website: www.elitefurnaces.com