

# Data Sheet

## Model: TSV16/25/450

1600 °C Maximum Temperature,  
Single Zone Vertical Tube Furnace



### INTRODUCTION

The TSV16 furnace is a vertical-mounted tube furnace with a separate control console, ideal for most general laboratory thermal processing applications.

The TSV16/25 Single Zone Tube furnace can accept a work tube with internal tube diameter of 25 mm. TSV16/25 range offers a choice of three heated zone lengths of 180mm, 450 mm & 610 mm.

TSV16/25/450 comes with a heated zone length of 450 mm.

### SPECIFICATIONS

**Maximum Temperature:** 1600 °C

**Maximum Continuous Temperature:** 1600 °C

| These models are heated by silicon carbide rods arranged equally around and parallel with the work tube for fast heat up and excellent temperature uniformity.

| This furnace design requires the use of a separate work tube of a grade suitable for the maximum temperature rating of 1600 °C.

| A work tube is not supplied as an integral part of the furnace and therefore needs to be purchased with the order as it is an essential accessory.

| Temperature Sensor: 'R' type Thermocouple.

| Energy efficient, high quality, low thermal mass insulation.

| High end Microprocessor PID controller.

| Vertical models are supplied as separate control console remote from furnace, separated by 2 meters of flexible interconnecting cable.



## TSV16/25/450

| **External Dimensions (mm):** H x W x D    **Net Wt.**  
(Nominal) TSV16/25/450

| Supply: 400V– 3 Phase – 50/60Hz

| Energy rating: 5.0 kW (Average)

33.0 kW (Peak)

### OPTIONAL FEATURES

| Work tubes of various materials, lengths and diameters for use in the furnace.

| The work tubes are available for containment of atmosphere or protection against process contaminants.

| A variety of triple flange gas tight end seals for work tubes to allow processing under vacuum or gas atmospheres.

| Multi segment, multi program storage controllers.

| Over temperature protection controller.

| Thermal Radiation Plugs.

## Elite Thermal Systems Limited

Elite Court, 6 Stuart Road, Market Harborough, Leicestershire LE16 9PQ, UK

Tel: +44 (0)1858 469834 | E-mail: [contact@elitefurnaces.com](mailto:contact@elitefurnaces.com) | Website: [www.elitefurnaces.com](http://www.elitefurnaces.com)