

Data Sheet

Model: TSV16/50/610

1600 °C Maximum Temperature,
Single Zone Vertical Tube Furnace



INTRODUCTION

The TSV16 furnace is a vertical-mounted tube furnace with a separate control console, ideal for most general laboratory thermal processing applications.

The TSV16/50 Single Zone Tube furnace can accept a work tube with internal tube diameter of 50 mm. TSV16/50 range offers a choice of three heated zone lengths of 180mm, 450 mm & 610 mm.

TSV16/50/610 comes with a heated zone length of 610 mm.

SPECIFICATIONS

Maximum Temperature: 1600 °C

Maximum Continuous Temperature: 1600 °C

| These models are heated by silicon carbide rods arranged equally around and parallel with the work tube for fast heat up and excellent temperature uniformity.

| This furnace design requires the use of a separate work tube of a grade suitable for the maximum temperature rating of 1600 °C.

| A work tube is not supplied as an integral part of the furnace and therefore needs to be purchased with the order as it is an essential accessory.

| Temperature Sensor: 'R' type Thermocouple.

| Energy efficient, high quality, low thermal mass insulation.

| High end Microprocessor PID controller.

| Vertical models are supplied as separate control console remote from furnace, separated by 2 meters of flexible interconnecting cable.



TSV16/50/610

| **External Dimensions (mm):** H x W x D **Net Wt.**
(Nominal) TSV16/50/610

| Supply: 400V– 3 Phase – 50/60Hz

| Energy rating: 7.0 kW (Average)
29.0 kW (Peak)

OPTIONAL FEATURES

| Work tubes of various materials, lengths and diameters for use in the furnace.

| The work tubes are available for containment of atmosphere or protection against process contaminants.

| A variety of triple flange gas tight end seals for work tubes to allow processing under vacuum or gas atmospheres.

| Multi segment, multi program storage controllers.

| Over temperature protection controller.

| Thermal Radiation Plugs.

Elite Thermal Systems Limited

Elite Court, 6 Stuart Road, Market Harborough, Leicestershire LE16 9PQ, UK
Tel: +44 (0)1858 469834 | E-mail: contact@elitefurnaces.com | Website: www.elitefurnaces.com